



Link for latest version of target value sheet:

<https://www.human.de/target-value-sheets/>

Target Values / Sollwerte / Valores Asignados/ Valeurs Cibles

	Method Methode Método Méthode	SI Unit SI Einheit Unidad SI Unité SI	Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs	Unit Einheit Unidad Unité	Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs
HUMAN reagent kits							
Level 1							
RHEUMATOID FACTORS							
HumaLyzer photometer	Immunoturbidimetry Latex enhanced	kIU/l	21.1	12.2 - 30.0	IU/ml	21.1	12.2 - 30.0
HumaStar analyzer		kIU/l	22.9	16.0 - 29.7	IU/ml	22.9	16.0 - 29.7
Level 2 *							
RHEUMATOID FACTORS							
HumaLyzer photometer	Immunoturbidimetry Latex enhanced	kIU/l	74.6	65.2 - 84.0	IU/ml	74.6	65.2 - 84.0
HumaStar analyzer		kIU/l	84.7	59.8 - 110	IU/ml	84.7	59.8 - 110

*For concentrations above the measuring range, a technically reasonable dilution has to be carried out.
Repeat the measurement and multiply the result by the dilution factor.